

# Nutrient Advantage

## Test Method References



SOIL CHEMICAL METHODS	Analysis equipment	Reference	Rayment & Lyons method
Analyses performed on soil dried at 40°C and ground to <2mm using a hammer mill and sieve			
Active Carbon (0.033M KMnO <sub>4</sub> )	FIA	1	6E1
Ammonium Nitrogen & Nitrate Nitrogen (1M KCl)	FIA	1	7C2b
Chloride (1:5 water)	FIA	1	5A2b
Boron (Hot CaCl <sub>2</sub> )	ICP-OES	1	12C2
Calcium Carbonate Percentage (% CaCO <sub>3</sub> equivalent)	Titration	1	19A1
“Fizz” test (carbonate qualitative field test)	Visual assessment	1	19D1
Electrical Conductivity (EC, 1:5 water)	Combination electrode	1	3A1
pH (1:5 water)	Combination electrode	1	4A1
pH (1:5 CaCl <sub>2</sub> )	Combination electrode	1	4B4
Exchangeable Aluminium (1M KCl)	ICP-OES	1	15G1
Exchangeable Cations: Calcium, Magnesium, Potassium, Sodium (1M Ammonium Acetate)	ICP-OES	1	15D3
Exchangeable Cations - Calcium, Magnesium, Potassium, Sodium (1M Ammonium Acetate, with soluble salt pre-treatment)	ICP-OES	1	15D1
Exchangeable Cations - Calcium, Magnesium, Potassium, Sodium (Alcoholic NH <sub>4</sub> Cl at pH 8.5, with soluble salt pre-treatment, Tucker)	ICP-OES	1	15C1
Exchangeable Cations - Calcium, Magnesium, Potassium, Sodium, Aluminium (BaCl <sub>2</sub> /NH <sub>4</sub> Cl, Gillman & Sumpter)	ICP-OES	1	15E1
Mehlich 3 extractable elements: P, Ca, Mg, Na, K, Cu, Fe, Mn, Zn, B, S, Al	ICP-OES	1	18F1
Molybdenum (hot CaCl <sub>2</sub> )	ICP-MS	1	12E1
Organic Carbon (Walkley & Black)	UV-VIS	1	6A1
Phosphorus Buffer Index, PBI (requires Phosphorus Colwell test)	ICP-OES	1	9I2b
pH Buffer (Liming requirement, Mehlich single buffer)	Combination electrode	1	16C1
Phosphorus (DGT)	FIA	2	
Phosphorus (BSES, H <sub>2</sub> SO <sub>4</sub> )	FIA	1	9G2
Phosphorus (Bray I)	FIA	1	9E2
Phosphorus (Bray II)	FIA	3	
Phosphorus (Olsen)	FIA	1	9C2b
Phosphorus (Colwell)	FIA	1	9B2
Potassium (Colwell)	Flame AAS	1	18A1
Potassium (Nitric)	ICP-OES	1	18C1
Potassium (Skene)	ICP-OES	4	
Potentially Mineralisable Nitrogen, PMN (Requires Moisture test. 10 day TAT)	FIA	1	7D2b
Saturated Paste/Extract: Electrical Conductivity	Combination electrode	1	14B1
Saturated Paste/Extract: Calcium, Magnesium, Sodium, Potassium, Sulphur, Phosphorus	ICP-OES	1	14F1
Silicon (BSES, H <sub>2</sub> SO <sub>4</sub> )	FIA	1	13D1
Silicon (CaCl <sub>2</sub> )	FIA	5	
Sulphur (KCl 40°C)	ICP-OES	1	10D1
Sulphur (MCP)	ICP-OES	1	10B3

# Nutrient Advantage

## Test Method References



### SOIL CHEMICAL METHODS

Analyses performed on soil dried at 40°C and ground to <2mm using a hammer mill and sieve

	Analysis equipment	Reference	Rayment & Lyons method
Total Nitrogen (Kjeldahl digestion)	FIA	1	7A1
Total Nitrogen (Combustion, Elementar Vario MAX Cube)	Combustion	1	7A5
Total Carbon (Combustion, Elementar Vario MAX Cube)	Combustion	1	6B2b
Total (Acid digestion) Phosphorus, Aluminium, Cadmium, Calcium, Chromium, Copper, Iron, Lead, Magnesium, Manganese, Nickel, Potassium, Sodium, Sulphur, Zinc	ICP-OES	1	17B1
Trace Elements: Copper, Iron, Manganese, Zinc (DTPA)	ICP-OES	1	12A1
Zinc (HCl)	ICP-OES	1	12D1
Cation Exchange Capacity (CEC)	Calculation	1	15J1
Calcium / Magnesium ratio	Calculation	1	15M1
Percentage Aluminium Saturation	Calculation	1	15O1
Exchangeable Sodium Percentage (ESP)	Calculation	1	15N1
Electrical Conductivity of Saturated Extract (ECSE)	Calculation	1	15J1

### SOIL PHYSICAL METHODS

Analyses performed on soil 'as received'.

	Analysis equipment	Reference
Dispersion & Slaking (Index 1 to 16, Normal agricultural assessment)	Visual assessment	6
Dispersion (Classes 1 to 8, Construction or agricultural assessment)	Visual assessment	7
Moisture (at 105°C)	Oven/Balance	8
%Gravel	2mm Sieve & Balance	
%Sand (2mm – 0.02mm), %Silt (0.02 – 0.002mm, %Clay (<0.002mm) & Texture (Australian Texture Triangle)	Soil Hydrometer	9
%Course Sand (2mm – 0.2mm), %Fine Sand (0.2mm – 0.02mm), %Silt (0.02 – 0.002mm, %Clay (<0.002mm) & Texture (Australian Texture Triangle)	Soil Hydrometer & Wet Sieve	9
Soil Colour (colour group via comparison to soil Munsell charts)	Visual assessment	10
Texture (Northcote field technique)	Subjective assessment	11

### PLANT TISSUE, GRAIN, MANURE, COMPOST, MILL MUD/ASH

Analyses performed on plant material dried at 70°C and ground to <1mm

	Analysis equipment	Reference
Ammonium Nitrogen (1:125 water extract)	FIA	12
Nitrate Nitrogen (1:125 water extract)	FIA	13
Chloride (1:125 water extract)	FIA	14
Moisture (at 70°C)	Oven/Balance	15
Total Elements: Boron, Calcium, Copper, Iron, Magnesium, Manganese, Phosphorus, Potassium, Sodium, Sulphur, Zinc ( $\text{HNO}_3$ and $\text{H}_2\text{O}_2$ digestion)	ICP-OES	16
Total Elements: Cobalt, Molybdenum, Selenium ( $\text{HNO}_3$ and $\text{H}_2\text{O}_2$ digestion)	ICP-MS	16
Total Nitrogen (Kjeldahl digestion)	FIA	17
Total Nitrogen (Combustion, Elementar Vario MAX Cube)	Combustion	17
Total Carbon (Combustion, Elementar Vario MAX Cube)	Combustion	18

# Nutrient Advantage

## Test Method References



WATER METHODS	Analysis equipment	Reference	APHA methods
Aluminium	ICP-OES	19	3120 B; 3030 B
Bicarbonate Alkalinity	Titration	19	2320 B
Carbonate Alkalinity	Titration	19	2320 B
Total Alkalinity	Titration	19	2320 B
pH	Combination electrode	19	4500-H+B
Electrical Conductivity	Combination electrode	19	2510 B
Chloride	FIA	19	4500-Cl-G
Ammonium Nitrogen	FIA	19	4500-NH3-H
Nitrate Nitrogen	FIA	19	4500-NO3-I
Sodium	ICP-OES	19	3120 B; 3030 B
Potassium	ICP-OES	19	3120 B; 3030 B
Calcium	ICP-OES	19	3120 B; 3030 B
Magnesium	ICP-OES	19	3120 B; 3030 B
Sulfur	ICP-OES	19	3120 B; 3030 B
Iron	ICP-OES	19	3120 B; 3030 B
Boron	ICP-OES	19	3120 B; 3030 B
Copper	ICP-OES	19	3120 B; 3030 B
Zinc	ICP-OES	19	3120 B; 3030 B
Manganese	ICP-OES	19	3120 B; 3030 B
Phosphorus	ICP-OES	19	3120 B; 3030 B

## REFERENCES

1. Rayment, G.E. and Lyons, D.J. 2011). Soil Chemical Methods - Australasia, CSIRO Publishing. Melbourne, Australia.
2. Mason, S., McNeill, A., McLaughlin, M.J., and Zhang, H. 2010) Prediction of wheat response to an application of phosphorus under field conditions using diffusive gradients in thin-films (DGT) and extraction methods. Plant and Soil 337: 243–258.
3. Bray, R.H. and Kurtz, L.T. 1945, Determination of total, organic, and available forms of phosphorus in soils, Soil Science 59: 39-45.
4. Skene, J. K. M. 1967), 'Soil analysis as an aid to diagnosing deficiencies of phosphorus and potassium'. In "Proceedings Australian Plant Nutrition Conference, Melbourne". pp 146 – 53
5. Haysom, M. B. C. and L. S. Chapman. 1975. Some aspects of the calcium silicate trials at Mackay. Proc. Qld. Soc. Sugar Cane Technol. 42:117-122
6. Loveday, J. and Pyle, J. 1973) The Emerson dispersion test and its relationship to hydraulic conductivity. Division of Soils Technical Paper No. 15. CSIRO, Melbourne, Australia.
7. Australian Standard AS 1289.3.8.1-2006. Determination of Emerson class number of a soil. In 'Methods of testing soils for engineering purposes. Method 3.8.1: Soil classification tests - Dispersion'. Standards Australia.
8. Rayment, G.E. and Lyons, D.J. 2011). 2A1 Air-dry moisture content. In Soil Chemical Methods - Australasia, CSIRO Publishing. Melbourne, Australia.

# Nutrient Advantage

## Test Method References



9. Gee, G.W. & Or, D. 2002) Particle-size Analysis p255-293. In J.H. Dane & G.C. Topp ed) Methods of Soil Analysis – Physical Methods. Part 4. SSSA, Madison, WI, USA.
10. Munsell Soil Color Charts 1992). Macbeth, New York.
11. Northcote, K.H. 1984). A factual key for the recognition of Australian soils. Rellim Technical Publications, Adelaide, SA. 4th Edition. pp:26-30
12. Kalra, Y. P. 1998). Determination of Ammonium-Nitrogen in Plant Tissue. In Handbook of Reference Methods for Plant Analysis. CRC Press, Boca Raton, USA. pp 89 – 91.
13. Kalra, Y. P. 1998). Extractable Nitrate in Plant Tissue. In Handbook of Reference Methods for Plant Analysis. CRC Press, Boca Raton, USA. pp 85 – 88.
14. Kalra, Y. P. 1998). Determination of Chloride in Plant Tissue. In Handbook of Reference Methods for Plant Analysis. CRC Press, Boca Raton, USA. pp 111– 113.
15. Kalra, Y. P. 1998). Determination of Dry Matter Content of Plant Tissue. In Handbook of Reference Methods for Plant Analysis. CRC Press, Boca Raton, USA. pp 51 - 52.
16. G. M. Bryson et. al. (2014). Wet Acid Digestion Procedure Using Nitric Acid and Peroxide, in Plant Analysis Handbook III. Micro-Macro Publishing. Georgia, U.S.A. p 227.
17. Kalra, Y. P. 1998). Determination of Total Nitrogen in Plant Tissue. In Handbook of Reference Methods for Plant Analysis. CRC Press, Boca Raton, USA. pp 75 – 83.
18. Rayment, G.E. and Lyons, D.J. 2011). 6B2b Carbon. in ‘Soil Chemical Methods – Australasia’. CSIRO Publishing, Melbourne, Australia. p 65 – 92
19. American Public Health Association (APHA). 1995. Standard Methods for the Examination of Water and Wastewater. 19th Edition, APHA, Washington, DC.